SISTEMA DE TROCAR DESECHABLE

ENDOSPHERE

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ENDOSPHERE

INSTRUCTIONS FOR USE

PLEASE READ THE FOLLOWING INFORMATION CAREFULLY

Failure to follow the safety instructions may result in serious complications during surgery.
IMPORTANT : these instructions are intended to provide information on the use of the "Endosphere"
range of trocar products and are not to be used as a point of reference for surgical procedures.

	Read instruction for use		Does not contain Latex
	Ethylene Oxide Sterilization		Single Use
	Expiration Date :YYYY-MM		Do not re-sterilize
	Manufacturing Date :YYYY-MM		Rx ONLY: Federal law (USA) restricts this device to sale by or on the order of a physician.

INDICATIONS

The Endosphere trocar is prescribed for laparoscopic applications as a port of entry for instruments.

PERFORMANCE

The Endosphere trocar ensures a full seal thereby maintaining the pneumopentoneum, including during instrument insertion.

CONTENTS

The Endosphere surgical trocar is available in three diameters. It consists either of a blade trocar with shielded function, either of a pyramidal tip trocar or of a blunt-tipped trocar for "open coelo" procedures.The sleeve can be smooth or threaded, with or without insufflation valve.

RANGE

Trocar size	Seal valve type	For use with the following instrument sizes
5 mm	Single valve	4.0 mm to 5.7 mm
5 / 11 mm	Detachable, double valve universal reducer	5.0 mm to 11.7 mm
5 / 12 mm	Detachable, double valve universal reducer	5.0 mm to 13.0 mm

INSTRUCTIONS

- Before opening, make sure that the sterile packaging has not been damaged.
- The trocar has been packed with the insufflation valve in the open position.
- Remove the trocar and the sleeve from the shell, and then re-assemble by sliding the trocar into the sleeve.
- Caution: the trocar should be disarmed before it is inserted into the sleeve to prevent any risk of damage to the membrane.***
- Check that the trocar has been assembled correctly.

Inserting the pyramidal tip trocar:

- Make a sufficient incision in accordance with trocar diameter.
- Hold the trocar firmly,keeping the cannula and the pyramidal tip obturator together.
- Apply a constant downward pressure to insert the obturator and cannula into the incision under visual control.
- When the cannula is in place in the wall, the pyramidal tip obturator is removed, and if necessary, it can be used with the second cannula available in the set.

Inserting the dual-safety blade trocar:

- Arm the trocar by pressing the colour-coded push button on the top of the Endosphere trocar.
- Apply a constant downward pressure to insert the trocar and sleeve into the incision.
- Once the abdominal wall has been penetrated, the shielded trocar function enables the cutting edge to be withdrawn while the shield is moving into place.
- Maintaining a steady pressure on the push button applies a scalpel function to the trocar.
- Caution: this operation cancels the safety mechanism feature.***
- Withdraw the trocar.

Inserting the Endosphere for open laparoscopy:

- Make a sufficiently long umbilical incision in accordance with open laparoscopy techniques.
- Place the trocar and anchoring device over the incision.Next,use the two wing screws to tension the anchoring device and adjust the sleeve height.
- Holding the membrane and the trocar head together,insert the blunt-tip obturator into the umbilical incision,down to the peritoneal cavity, according to the open laparoscopic technique.
- The assembly is secured by two sutures on the anchoring device's suture stays.
- Withdraw the blunt-tipped trocar.

SAFETY INSTRUCTIONS

- "Federal (USA) law restricts this device to sale by or on the order of a physician".**
This instrument is intendeded for use by practitioners with suitable training in laparoscopic surgery.
- Endosphere trocar sleeves must only be used with Endosphere trocars.
Users must have a complete understanding of the risks involved in electro surgical procedures.
Consult medical papers if necessary.
 - All instruments used during the procedure should be checked for compatibility.
 - Instruments should be inserted along the trocar's centre line.Angled instruments should be handled with care to avoid wrenching the trocar valve free.

STORAGE

Keep the sterile packing unopened in a cool, dry and well-ventilated place away from extreme temperatures and humidity.

LIFE SPAN

The pack contents will remain sterile as long as the package has not been tampered with and are to be used before the use-by date shown on the package. Do not use if the package has been opened or damaged.

SINGLE-USE TROCAR: DO NOT RESTERILIZE

The manufacturer accepts no responsibility for problems resulting from the inappropriate use of products. Prices, details and availability may be changed without prior notice.

In the event of reuse or re-sterilization of these instruments designated for single use, the retention of physical properties cannot be guaranteed. Re-sterilization could cause the deformation of plastic parts thereby obstructing instruments being introduced, preventing the circulation of insufflationary gas, and leading to a weakening of the seal.

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